SABLE

INDUSTRIES INC.

ISO 13485:2003 CERTIFIED • HEALTH CANADA REGISTERED

Instructions

Name: Sable Access Prophy Clean
Part Number: 1600015

Sable – Hygienist Handpiece (Midwest RDH Style)

Technical Facts


ISO U-type Doriot Motor Maximum rpm: 0- 5000 min
Attachment: accepts all disposable doriot prophy angles
Air Requirements: Clean filtered moisture free air with recommended pressure of 30- 45 psi.
Do not exceed 60 psi.

Operation

Connect and disconnect prophy angle from motor:

Connect: Press the angle firmly into the motor as shown in the picture.
The motor must fully engage and fit tightly with attachment.

Disconnect: The motor must not be moving before disconnecting the prophy angle.
Pull the motor and angle apart with a firm pull.

Cleaning and Care

Your Sable hygienist handpiece is a high quality precision instrument. Incorrect maintenance and care can shorten the life of this product. We recommend for cleaning and lubrication FDA approved Spray or DCL-90.

External Cleaning

Remove angle from attachment and disconnect from hose and clean the external surface thoroughly with a sponge or gauze using warm tap water.

DO NOT IMMERSE INSTRUMENT. Wipe dry with a clean cloth.

Clean/Lubricate Internal Surface

Lubricate before every sterilization. Use a combination FDA approved cleaner/lubricant and spray into the drive air tube. Also put lubricant inside of chuck where the prophy angle is inserted. Run motor on air source for 30 seconds so that all excess lubricant is expelled. This ensures all internal parts of motor are completely lubricated.

2 Hole

4 Hole

Drive air Hole

Weekly – Place a few drops of lubricant in the swivel junction in the middle of the handpiece and twist swivel to work lubricant into junction and wipe clean.

AUTOCLAVING PROCEDURE:

The handpiece does not come sterilized from the manufacturer. It is also recommended that the handpiece or motor be sterilized after each use. FDA/Health Canada autoclaving or sterilization packaging is recommended to preserve item in sterile condition.

1. Wipe off debris on the entire tool surface with alcohol-soaked cloth. Phenolic Alcohol, up to 3% glutaraldehyde based on the agents containing acetone or chlorine.
2. Lubricate (See Lubrication)
3. Follow the instructions for your particular Autoclave or Sterilizer

CAUTION

DON'T use with burs or metal prophy angles.
DON'T immerse instrument in water or use ultrasonic cleaner.
DON'T use any type of disinfectant, chemical or soap on instrument.
Use of chlorine products, aldehydes, etc. will damage handpiece and void all warranties.
DON'T operate motor while inserting/removing attachments
DON'T exceed 275°F (29 psi) / 135°C (2 bar) in autoclave.
DON'T dry heat or heat transfer sterilize

Sable Industries Inc.
100 Campbell Ave., Unit 5 Kitchener ON
Canada Phone 519-579-9323 Fax 519-579-9324
www.sableindustriesinc.com

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