**MAINTENANCE & CLEANING MANUAL**

1600111  
SABLE 20K E-Type 4-Hole Mini Air Motor

1600110  
SABLE 20K E-Type 4-Hole Air Motor

1610110  
SABLE 20K E-Type 4-Hole Air Motor with Internal Air/Water

1600102  
SABLE 5K E-Type 4-Hole Air Motor

1600005  
SABLE 5K U-Type Hygienist  
Prophy Swivel Handpiece

**SPECIFICATIONS**

| Model                | 1600111<br>SABLE 20K<br>E-Type 4-Hole<br>Mini Air Motor | 1600110<br>SABLE 20K<br>E-Type 4-Hole<br>Air Motor | 1610110<br>SABLE 20K<br>E-Type 4-Hole<br>Air Motor with<br>Internal Air/<br>Water | 1600102<br>SABLE 5K<br>E-Type 4-Hole<br>Air Motor | 1600005<br>SABLE<br>5K U-Type<br>Hygienist<br>Prophy Swivel<br>Handpiece |
|----------------------|---|--|---|---|--|
| Handpiece Type       | Slow Speed Air Driven Motor                             | Slow Speed Air Driven Motor                        | Slow Speed Air Driven Motor   | Slow Speed Air Driven Motor                       | Slow Speed Air Driven Motor  |
| Attachment Type      | ISO 7785-2 Type and OEM type                            | ISO 7785-2 Type and OEM type                       | ISO 7785-2 Type and OEM type  | ISO 7785-2 Type and OEM type                      | U-Type Heads Only  |
| Rotation Speed       | Up to 20,000 RPM  | Up to 20,000 RPM                                   | Up to 20,000 RPM  | Up to 5000 RPM                                    | Up to 5000 RPM   |
| Air Pressure         | 3 bar =43.5 PSI   | 3 bar =43.5 PSI                                    | 3 bar =43.5 PSI   | 3 bar =43.5 PSI                                   | 3 bar =43.5 PSI  |
| Weight               | 90 grams  | 90 grams   | 85 grams  | 110 grams   | 75 grams   |
| Sterilization        | Autoclave 135°C @3min, dry time =16min                  | Autoclave 134°C @8min, dry time=20 min             | Autoclave 134°C @8min, dry time=20 min  | Autoclave 132°C @10min, dry time=40 min           | Autoclave 132°C @10min, dry time=40 min                                  |
| Sterilization Cycles | 250 + as per ISO 7785-1                                 | 250 + as per ISO 7785-1                            | 250 + as per ISO 7785-1   | 250 + as per ISO 7785-1                           | 250 + as per ISO 7785-1  |

**! SPECIAL PRECAUTIONS:**

- Always exercise extreme caution when operating your Motor or Prophy Handpiece.
- Always inspect your Motor or Prophy Handpiece thoroughly before each use.
- Always use PPE (Proper Protective Equipment) when operating and cleaning your Sable Motor or Prophy Handpiece.
- See your Occupational Health and Safety Act or other safety ordinances for specific instructions.
- This is a precision instrument. Do not drop or apply excessive force when handling.
- If product is dropped your dealer should be contacted, it is recommended that you do not use.
- Contact your dealer right away if any problem should occur with any of your Sable products.
- Sterilize before first time use and before every subsequent use thereafter.
- Indication for use: Used by a trained dental practitioner in conjunction with other attachments for the cleaning of teeth.

**HOW TO USE:**

- Attachment of E-Type Motor or Prophy Handpiece. Align the pins on the Motor or Prophy Handpiece with the holes in the Air tubing turn the collar on the tubing clockwise until snug.
- To remove the Motor or Prophy Handpiece, turn the collar on the tubing counter clockwise until loose.

**CYCLES:**

The Motor and Prophy Handpiece can be subject to 250 cycles = (Performed Sterilization Process) per EN ISO 5.10 before showing signs of deterioration.

**CLEANING: Preparation and Cleaning before Sterilization:**

Clean and lubricate the handpiece as soon as possible after treatment. Failure to properly maintain the handpiece may cause infection, overheating leading to burn injuries.

**! DO NOT IMMERSE INSTRUMENTS IN DISINFECTANT SOLUTION OR CLEAN IN ULTRASONIC DEVICES. MANUAL CLEANING ONLY.**

**1. PREPARATION:**

- Always wear protective gloves, a mask, and protective goggles for safety purposes and to minimize the risk of infection.
- Wipe the exterior of the handpiece/accessory using a cleaning wipe.
- Operate the handpiece at the chair-side for at least 20 seconds to purge fluids in the handpiece (not applicable for air motors).
- Remove the handpiece, carry the handpiece to the decontamination area and proceed to step 2.

**2. CLEANING AND DRYING:**

*\*Wash the handpiece and irrigation nozzle under running water following the procedures below. (Water condition:<38°C (100°F), >3.5L/min; water should be same quality as drinking water)*

- Clean the external surfaces of the handpiece/accessory with a soft bristled toothbrush for 15 seconds or more.
- Under appropriate lighting (natural office lighting), inspect the handpiece/ accessory for blood or other debris. If any visible debris remain, repeat the process until it is visually clean.
- Wipe the exterior of the handpiece/accessory clean using a cloth moistened with disinfectant or ethanol with a concentration of 60-90% (ie. Isopropyl).

**LUBRICATION:**

Recommended to use Sable Branded autoclavable handpiece lubricant.

- Remove Motor or Prophy Handpiece from airline.
- Using spray lubrication can, oil dropper or lubrication station, apply lubricant in the drive air tube (smaller of the two tubes at the rear of the motor).
- Reattach Motor or Prophy Handpiece to airline if not using a lubrication station.
- Run entire tool at normal speed for 10 seconds.
- Wipe off excess oil from entire tool surface.

**! AUTOCLAVING/STERILIZATION**

**IMPORTANT NOTE:** This handpiece or motor do not come sterilized from the manufacturer. It is also recommended that the handpiece or motor be sterilized after each use. Sable Branded autoclaving or sterilization packaging is recommended to preserve item in sterile condition. Please note that after sterilization, handpiece should be used at your earliest convenience to make certain that sterilization is maintained.

- Clean before sterilization (see CLEANING above)
- Lubricate before sterilization (see LUBRICATION above).
- For part number **1600111** set parameters for Moist Heat Sterilizer. PreVac cycle 135 degrees Celsius for 3 minutes and 16 minute dry time.
- For part numbers **1600110, 1610110** set parameters for Moist Heat Sterilizer. PreVac cycle 134 degrees Celsius for 8 minutes and 20 minute dry time
- For part numbers **1600005, 1600102** set parameters for Moist Heat Sterilizer. PreVac cycle 132 degrees Celsius for 10 minutes and 40 minute dry time.

**WARRANTY INFORMATION:**


Sable Slow Speed Motor or Prophy Handpiece are warranted to be free of manufacturing defects for a period of 12 months\* from date of purchase for the standard Motor or Prophy Handpiece.

*\*Warranty is conditional on proper use and maintenance of the Motor or Prophy Handpiece. Warranty is limited to repair or replacement at manufacturers' discretion. Under no circumstances will the manufacturer be held responsible for incremental or consequential damages. Please contact point of purchase authorized distributor for warranty repair.*

**DISCLAIMER:**

*Sable Industries Inc. recommends that the sterilization process outlined in the Directions for use (DFU) described in this document be followed precisely. Sable Industries Inc. is not responsible for any products supplied in which the end user does not follow the DFU.*

DEALER

Manufactured by  
SABLE INDUSTRIES INC.  
100 Campbell Ave. Suite 5  
Kitchener, ON N2H 4X8, Canada  
TEL: 1-800-368-8106  
FAX: 1-519-579-9323  
Email: info@sableindustriesinc.com  
www.sableindustriesinc.com

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1600111  
SABLE 20K E-Type 4-Hole Mini Air Motor

1600110  
SABLE 20K E-Type 4-Hole Air Motor

1610110  
SABLE 20K E-Type 4-Hole Air Motor  
with Internal Air/Water

1600102  
SABLE 5K E-Type 4-Hole Air Motor

1600005  
SABLE 5K U-Type Hygienist  
Prophy Swivel Handpiece

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