# SABLE INDUSTRIES INC.✿

## SABLE SLOW SPEED HEADS MAINTENANCE & CLEANING MANUAL

**5,000 RPM** - 1600401, 1600402, 1600403, 1600405, 1600406, 1600408, 1600409, 1600413, 1600501, 1600502, 1600503, 1610820, 1610821

**20,000 RPM** - 1600404, 1600407, 1600410, 1600412, 1600421, 1600504, 1600507 1610810

### **SPECIAL PRECAUTIONS:**

1. Always exercise extreme caution when operating your Slowspeed.
2. Always inspect your Slowspeed thoroughly before each use.
3. Always use PPE (Proper Protective Equipment) when operating and cleaning your able Slowspeed.
4. See your Occupational Health and Safety Act or other safety ordinances for specific instructions.
5. This is a precision instrument. Do not drop or apply excessive force when handling.
6. If product is dropped your dealer should be contacted, it is recommended that you do not use.
7. Contact your dealer right away if any problem should occur with any of your Sable products.
8. Sterilize before first time use and before every subsequent use thereafter.
9. Indication: in conjunction with other attachments for the cleaning of teeth.

### **HOW TO USE:**

**Caution:** For maximum speed, selection ring on motor must be fully in F or R position.

1. Attachment of Head to E-Type Contra Angles. Insert the head and shaft (in the head) in the crown of the contra angle. Ensure the crowns are aligned so the teeth on the contra angle and the head mesh together. Twist (turn) collar on head to the right to tighten the contra angle and the head together, (only finger tight).
2. Removal of Head from the Contra Angles. Twist (turn) collar on head to the left to loosen. Ensure the head transmission shaft is taken out and placed back into head.

### **CYCLES:**

The Attachment/Head can be subject to 250 cycles = (Performed Sterilization Process) per EN ISO 14457 requirement 5.10 before showing signs of deterioration.

### **FEATURES & SPECIFICATIONS:**

1. Head conforms to ISO 14457 specifications.
2. Accepts all E-Type (ISO 3964 compatible) attachments.
3. Connection with Attachments: ISO Standard.
4. Method of Sterilization: Autoclave.

### **CLEANING: Preparation and Cleaning before Sterilization:**

*\*Clean and lubricate the handpiece as soon as possible after treatment. Failure to properly maintain the handpiece may cause infection, overheating leading to burn injuries.  DO NOT IMMERSE INSTRUMENTS IN DISINFECTANT SOLUTION OR CLEAN IN ULTRASONIC DEVICES*

#### **1. PREPARATION:**

- a) Always wear protective gloves, a mask, and protective goggles for safety purposes and to minimize the risk of infection.
- b) Remove the bur. (if applicable)
- c) Wipe the exterior of the handpiece/accessory using a cleaning wipe.
- d) Operate the handpiece at the chair-side for at least 20 seconds to purge fluids in the handpiece. (not applicable for air motors)
- e) Remove the handpiece, carry the handpiece to the decontamination area and proceed to step 2.

#### **2. CLEANING AND DRYING:**

*\*Wash the handpiece and irrigation nozzle under running water following the procedures below. (Water condition: <38°C (100°F), >3.5L/min; water should be same quality as drinking water)*

- a) Clean the external surfaces of the handpiece/accessory with a soft bristled toothbrush for 15 seconds or more.
- b) Clean around the bur insertion hole for 15 seconds or more, using an interdental brush (wire thickness 0.7mm) that conforms to size 4 in ISO 16409. (if applicable)
- c) Under appropriate lighting (natural office lighting), inspect the handpiece/accessory for blood or other debris. If any visible debris remain, repeat the process until it is visually clean.
- d) Wipe the exterior of the handpiece clean using a cloth moistened with disinfectant or ethanol with a concentration of 60-90% (ie. Isopropyl)

### **LUBRICATION:**

#### **USE ONLY AUTOCLAVABLE HANDPIECE LUBRICANT.**

1. Leave head attached to contra angle; remove from motor
2. Using a spray lubrication can or lubrication station, apply lubricant into contra angle.
3. Reattach entire contra angle and head to air motor.
4. Run entire tool at normal speed for 10 seconds.
5. Wipe off excess oil from entire tool surface.

### **AUTOCLAVING PROCEDURE:**

Sterilize in accordance with ISO17665-1:2006. Suitable autoclaving packaging is recommended to preserve item in sterile condition. Once heads has been sterilized it should be used as at earliest convenience to ensure sterilization is maintained.

1. Clean before sterilization (See Cleaning).
2. Lubricate before sterilization (See Lubrication).
3. For part numbers, **1600401, 1600402, 1600403, 1610404, 1600405, 1600406, 1600407, 1600408, 1600409, 1600413, 1600501, 1600502, 1600503, 1600504, 1600507** set parameters for Moist Heat Sterilizer. PreVac cycle 132 degrees Celsius for 10 minutes and 40 minute dry time.
4. For part numbers **1600410, 1600412, 160421, 1610810, 1610820, 1610821**, sterilization parameters for moist heat sterilizer pre vac cycle to 134 degrees Celsius for 8 minutes dry time 20 minutes.
5. Use Sterilized head as soon as possible to ensure sterilization is maintained

### **WARRANTY INFORMATION:**

The Sable Slow Speed Head is warranted to be free of manufacturing defects for a period of 3 months\* from date of purchase for the standard heads and 12 months on the Premium Heads.

*\*Warranty is conditional on proper use and maintenance of the Head. Warranty is limited to repair or replacement at manufacturers' discretion. Under no circumstances will the manufacturer be held responsible for incremental or consequential damages. Please contact point of purchase authorized distributor for warranty repair.*

**Disclaimer:** Sable Industries Inc. recommends that the sterilization process outlined in the Directions for use (DFU) described in this document be followed precisely. Sable Industries Inc. is not responsible for any products supplied in which the end user does not follow the DFU.

**Sable Industries Inc.**

100 Campbell Ave. Suite 5 Kitchener, ON N2H 4X8, Canada

See Sables Website for other Sable products and product information

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