
THE AUTOSCALER

ULTRASONIC DENTAL SCALER

DIRECTIONS FOR USE (DFU)

**PLEASE READ THOROUGHLY BEFORE OPERATING THIS DEVICE
THE INTENDED USE OF THE AUTOSCALER DEVICE IS
FOR DENTAL PROPHYLAXIS TREATMENT.**

**THIS DEVICE SHOULD BE USED ONLY BY LICENSED
DENTAL PROFESSIONALS TRAINED IN THE
PROPER USE OF ULTRASONIC SCALING DEVICES**

FEATURES

THE AUTOSCALER is designed for *Dental Prophylaxis Treatment* when used with the appropriate ultrasonic scaling inserts. If you have purchased **THE AUTOSCALER** 30,000Hz model you ***must*** use the **30K frequency insert**. If you have purchased **THE AUTOSCALER** 25,000Hz model you ***must*** use the **25K frequency insert**.

Inserts of different frequency ratings are not interchangeable.

The Autoscaler accepts all major brands of scaling tips.

THE AUTOSCALER name is owned by Sable Industries Inc. of Canada

THE AUTOSCALER is acceptable for use in scaling and cleaning surfaces of teeth in conjunction with other prophylactic measures. Revision: March 16, 2007

How to contact us:

Sable Industries Inc.

100 Campbell Ave. Suite 5

Kitchener Ontario

Canada N2H 4X8

519-579-9323 or toll free 800-368-8106

Recommendations by the manufacturer contained within this instruction booklet are recommendations only. Other options may be available or mandated by the professional or regulating agencies.

RESPONSIBILITIES

It is the responsibility of the dental professional using this device to be fully trained in the proper use of ultrasonic scalers. The operator must be observant of possible changes in performance during a procedure and immediately discontinue use of this product, and determine the cause of such changes in performance before continuing dental procedures. Performance varies slightly between manufacturers of ultrasonic scalers, so the operator must become familiar with the operating parameters, and adjustments of this particular device. Performance varies between Ultrasonic Inserts due to the style Insert, model Insert, and from the age and/or modifications of the insert.

SPECIAL PRECAUTIONS

- ✓ **NOTE:** The flexible **Water Input Hose** is subject to deterioration by aging, chemicals in the water supply, excessive flexing, and pressure fluctuations. As an added precaution, the water supply should be turned **OFF** at the end of the day. **REMEMBER** to turn the water back **ON** prior to use the next day.
- ✓ **CAUTION: THE AUTOSCALER** should **NOT** be used for restorative dental procedures involving the **CONDENSATION OF AMALGAM**.
- ✓ **CAUTION: THE AUTOSCALER** should **NOT** be used on or near anyone fitted with a **CARDIAC PACE MAKER**.
- ✓ **CAUTION: DO NOT** remove the **COVER** from the unit. There are **NO** user serviceable parts inside. Allow only qualified technicians to service **THE AUTOSCALER**.
- ✓ **CAUTION: DO NOT AUTOCLAVE** the detachable **HANDPIECE** assembly. A Federally approved non-immersion **WATER BASED** cold sterilization solutions **MAY** be used on the detachable handpiece assembly however **DO NOT** get sterilization solutions on the electrical contacts of the Handpiece **CONNECTOR** as this will cause oxidation, and a poor electrical connection. **DRY ALL SURFACES** with clean Air or clean paper towels.

STERILIZATION

- ✓ Manufacturer recommends that all operators of this ultrasonic unit wear at least one pair of **SURGICAL RUBBER GLOVES**, and for extra prevention wear double pairs of gloves to isolate the operator from infection, and cross contamination of other patients.
- ✓ It is recommended that complete **EYE** and **FACE** protection be used during the use of this unit to prevent Saliva Spray from being inhaled, and from entering the operators eyes.
- ✓ Operator should be careful not to **PUNCTURE THE SKIN** with any insert tip containing saliva from patients. Insert tips should be sterilized after each patient according to the ***insert manufacturers instructions.***
- ✓ Manufacturer recommends **DISPOSABLE HANDPIECE SHEATHES** that keep contamination from getting on the handpiece assembly. Dispose sheathe after each patient. Ask your local dealer for **HIGH SPEED LONG** sheathes.

INTRODUCTION


THE AUTOSCALER Ultrasonic Scaler is a precision-engineered instrument, constructed with only the finest of materials, and requiring no special maintenance. By taking the following precautions your unit will give you many years of service.

1. Do not place unit on or near any source of heat.
2. Clean the outer surfaces of **THE AUTOSCALER** with a Federally approved non-immersion WATER BASED cold sterilization solution. Dry immediately with clean air or clean paper towels.
3. Unit should be located where a normal amount of air will circulate freely around the case when in use.
4. The handpiece, and Quick Connector are very strong but care should be taken **NOT TO PULL** with excessive force.

(DESCRIPTION): Basically, the unit is broken down into **THREE** major components, ***the electronic generator, the water system, and the handpiece assembly.***

ELECTRONIC GENERATOR AND WATER SYSTEM

The electronic generator produces the frequency, and the power to activate a wide variety of interchangeable **Sable Scaling TIPS**. The water system provides precise control over the water (lavage) flow, which serves to wash debris from the work area, and to cool the laminations of the insert tip. The following controls are found on the front panel of **THE AUTOSCALER**.

- ✓ **TOGGLE SWITCH:** To turn unit On and Off.
- ✓ **INDICATOR L.E.D. LAMP:** (above the switch) which lights when the Toggle Switch is set to the **ON** position, indicating the unit is **ready for operation**.
- ✓ **POWER CONTROL KNOB:** NOTE: FOLLOW INSERT MANUFACTURERS COMPLETE INSTRUCTIONS, Linear Power Control knob adjusts the output power to the insert tip. The **BLUE** shaded area is for use when using modified **Periodontal inserts**. It is recommended the Power **not** be set higher than **12:00** even if no mist appears at the tip.
- ✓ **WATER / LAVAGE CONTROL KNOB:** For adjusting the flow of water to the handpiece, and insert tip. To obtain the **NECESSARY** water flow of at least **(18cc - 30cc+/minute)** - It may require **ONE or MORE** complete rotations of the knob in either direction. **HIGHER** water flow to the **RIGHT**, and **LOWER** water flow to the **LEFT**.
- ✓ **HANDPIECE PANEL RECEPTACLE:** Accepts the **MALE** side of the *Handpiece connector with lock ring*.

HANDPIECE ASSEMBLY


The handpiece assembly is used to accept all brands of insert tips. It is connected to **THE AUTOSCALER** front panel by a Quick Connector. The **open end of the handpiece will accept all brands of inserts, Inserts are frequency specific.**

- ✓ 30K tips measure (Shorter style approx. 6.25 inches long).
- ✓ 25K tips measure (Longer style approx. 7.25 inches long).

(25K) inserts will NOT work in your 30K unit.

(30K) inserts will NOT work in your 25K unit.

As illustrated both frequency inserts look identical except for length.

(NOTE): Inserts are extremely **FRAGILE**. If one is dropped, bent, or scratched it **WILL NOT** cavitate (*i.e. will non function*) or function by operating at a higher than normal temperature, and should be replaced **immediately**. **(NOTE):** Insert tips manufactured by other manufacturers may not achieve optimal results, and could run at a **higher temperature or not tune efficiently** in **THE AUTOSCALER**. **Insert tips should be routinely replaced according to insert manufacturers recommendations.**

ELECTRICAL REQUIREMENTS


A **grounded** wall receptacle is required supplying a voltage range of **110-125 volts 50/60 Hz A.C.**, or (*optional voltage range of*) **220-240 volts 50/60 Hz A.C.** If a 3 to 2 prong adapter is used, the ground wire **MUST** be connected to the base plate screw. **Have an electrician test your outlet** for correct wiring polarity, and a safe ground. Do not assume the ground is complete. **If you plug THE AUTOSCALER in a GFI (Ground Fault Interrupt) outlet it may trip the socket due to the high sensitivity of the socket. If this occurs plug THE AUTOSCALER into a non-GFI socket.**

INSTALLING THE HANDPIECE ASSEMBLY

1. Align the electrical contacts of the **MALE** connector with the same on the **FEMALE** panel connector.
2. Gently insert the connector. It should slide in with very little effort. If the connector does **NOT** align, **DO NOT** force it. Turn the **CONNECTOR BODY** in either direction until the connector correctly aligns, and goes in.
3. Once the connector has been completely inserted, screw the **OUTSIDE LOCK RING** to the **RIGHT** with light finger pressure until it stops. **DO NOT** over tighten the lock ring. **Be careful not to cross thread the lock ring.**
4. To **REMOVE** the connector, unscrew the **LOCK RING** to the left, and gently pull outward grasping the entire **CONNECTOR BODY**.

THE AUTOSCALER INSTALLATION AND OPERATION

1. **(NOTE): DO NOT use THE AUTOSCALER until all of the following steps have been completed.**
2. Upon receipt of your unit carefully unpack it and see that the Toggle Switch is in the **OFF-O** position.
3. The Water Control Knob should be in the **HIGH** position. **(Note: It is normal for the Knob to turn past the High position several rotations which will increase water flow substantially.)**
4. Set the **POWER KNOB** all the way to the **LEFT** in the **LOW** position.
5. Plug the **MALE** Water Quick Connector into the Dental operator or have a plumber install a SHUT OFF style FEMALE connector. **A minimum of (20 P.S.I.) water input pressure**, and a maximum of (60 P.S.I.) is needed from the operator for proper operation of **THE AUTOSCALER**. **NOTE: It is recommended that water line filtration be used to remove larger particles before entering the scaler.**
6. Plug the **POWER CORD** into the wall (*see electrical requirements*). Press the toggle switch to the **ON-I** position. The **L.E.D.** light above the switch should be ON ("I" position).
7. If your device is equipped with an **AIR SWITCH ACTIVATOR** the air supply must be attached to the rear of the device, and a **minimum of (20 PSIG)** is necessary to activate the air switch. **DO NOT REMOVE THE AIR**

RESTRICTOR FITTING. The restrictor relieves the air input line if **accidentally charged** from another air source without a purge cycle. The scaler air line will bleed off to zero pressure protecting **THE AUTOSCALER.**

8. Hold the handpiece in a **VERTICAL** upright position, and step on the **foot switch**. Within **10-15 seconds** water should come out of the open end of the handpiece **IF NOT, QUICKLY turn WATER KNOB clockwise 1-2 rotations to get water flow.** Some **AIR** will escape at the same time. This step is **IMPORTANT** because trapped air will cause **HOT SPOTS (usually noted by a loud screeching sound)** which can affect the performance of **THE AUTOSCALER** or permanently damage the handpiece assembly.

9. **(NOTE):** Inserts and handpiece will **BURN OUT** in **45-60 seconds** in the **ABSENCE OF WATER.** A **(18cc - 30cc+/min)** flow is the normal flow range, and necessary **at all times.** **NEVER** touch the Cavitating/Vibrating end of the insert tip to any **FLESH** such as the cheek, tongue, or gums unless performing special periodontal procedures. This will cause **a friction burn** to the patient. Always keep the insert in motion. Do not rest the handle or insert on the lips during a procedure. The **REASON** for this is due to the **25,000 or 30,000 back and forth strokes of the insert every second.** This is the principal idea of ultrasonic cleaning activity. The insert tip should touch **PRIMARILY** the calculus on the tooth surface, and **NOT** the enamel of the tooth surface.
10. Select an **INSERT.** Moisten the rubber **O RING** before inserting into the handpiece. While holding the handpiece **UPRIGHT,** place the insert in the **OPEN END** of the handpiece, and with a gentle **PRESS and TURN** motion, snap the insert in fully. **DO NOT** force the insert.
11. **(NOTE):** Manufacturer recommends disposable handpiece sheathes be used over the handpiece assembly. Place the sheathe over the handpiece. While holding the sheathe in place, puncture the **INSERT** through the plastic at the open end of the handpiece, and snap the insert in fully according to **STEP F.** Dispose sheathe after each patient.

Page Number 8

12. **(NOTE):** Due to atmospheric exposure, and exposure to sterilization solutions, the **O RING** on the insert tip may become swollen, brittle, or crack, and require replacement.
13. While holding the handpiece upright, press the Foot Switch, and water should begin to exit the insert. Allow the **AIR** bubbles to escape for about **10-20 seconds**. Turn the **POWER KNOB** to the **RIGHT** in the Medium position, and the insert tip should Cavitate/Vibrate with a strong mist. Turn the **WATER KNOB** to the right if there is **No mist or Weak misting**. If **Modified Periodontal inserts** are used, set the **POWER** knob in the **BLUE** shaded area. This power level is designed to be used with these **fragile** specialty inserts. Inserts may give a drip and mist combination. This is normal for these inserts, and is actually preferred by most Periodontal professionals due to the water drip being ultrasonically activated within **0.5mm** of the tip assisting in killing plaque.
14. The application of the **POWER KNOB** is for controlling the Output Power to the insert tip in a smooth linear action. The operator can select the **EXACT** power needed for cleaning. Set power knob **LOW** to **MEDIUM** for cooler handpiece/insert temperature.
15. If your insert will not Cavitate/Vibrate, the problem could be the **INSERT**. Before requesting service **TRY ANOTHER INSERT**. **Make sure you have the CORRECT FREQUENCY insert**. If you purchased a **30K** 30,000 Hz model, a sticker on the front panel or rear panel will indicate that you can use only 30K inserts. **(30K inserts will not function in a 25K unit, and 25K inserts will not function in a 30K unit)**.
16. The **HIGH** power setting is for applications where an *excessive amount of calculus* is built up. The Low power setting can be used for periodontal *Subgingival curettage, and root planing* procedures. Most scaling should be done only on the **LOW to MEDIUM** power setting.

UNIT IS NOW READY FOR USE


SHUT DOWN PROCEDURE

1. **IMPORTANT:** Turn off the water supply to **THE AUTOSCALER** by **disconnecting** the **QUICK WATER CONNECTOR**, or **shutting off** the main operator water supply. **REMEMBER** to turn the water back on the next day.
2. Press the **Toggle Switch** to the **OFF-O** position. The light should be off.

WHAT TO DO IF UNIT STOPS WORKING

If the unit does not function, proceed with the following:

1. If the power L.E.D. light above the switch is not lit, replace the **FUSE** on the rear panel with a "**LITTLE FUSE**" **2 amp Fast Blow style**. If the fuse opens again **DO NOT** attempt to replace the fuse as this indicates a serious condition that must be examined by the factory. Call for technical service **(519) 579-9323** or ask your local dealer to look up our current phone number and address in their records, or the **Proofs Magazine** listing of dental manufacturers.
2. Check the plug. **Is it all the way in the receptacle?**
3. **Try another insert.** (Make sure it is the correct frequency **25K or 30K** depending on your model.) If you have a 30K model a **30K sticker** will be on the front panel in the lower right corner or marked on the back of the unit usually on the serial number plate. **The most common mistake when purchasing a new unit is selecting the wrong frequency Cavitron insert.**
4. If water flow is **insufficient**, check water pressure at the operator for minimum of **20 PSI**. The **Insert Water SLOT** could be clogged with debris. Try another insert. See if **WATER QUICK CONNECTOR** is **locked fully** in place.
5. If the handpiece handle or insert operates **HOT: REDUCE POWER LEVEL and INCREASE WATER FLOW.**
6. **NEVER** touch the Cavitating/Vibrating end of the insert tip to any **FLESH** such as the cheek, tongue, or gums unless performing special periodontal procedures. This will cause **a friction burn** to the patient. Always keep the insert in motion. Do not rest the handle or insert on the lips during a procedure. The **REASON** for this is due to the **25,000 or 30,000 back and forth strokes of the insert every second**. This is the principal idea of ultrasonic cleaning activity. The insert tip should touch **PRIMARILY** the calculus on the tooth surface, and **NOT** the enamel of the tooth surface.
7. If **THE AUTOSCALER stops functioning or is intermittent** when the handpiece cord is moved, this means a wire in the handpiece cord is broken and must be immediately repaired or replaced. **DO NOT** continue to press the foot switch as there is a possibility the internal circuit **safety fuse** may blow and require replacement by a qualified technician.

Page Number 10

8. **If handpiece service is needed:** Turn the unit **OFF**. Carefully unscrew the outer lock ring of the handpiece assembly and gently unplug the handpiece from the front panel. You can call your dealer for assistance or **(519) 579-9323** for assistance.
9. If a **humming or high pitch sound** is noted from the console during routine use, this is a **normal** circuit resonance which will change as the POWER control knob is adjusted up and down, and will not affect the performance of **THE AUTOSCALER**.
10. A note about **Slimline™** and **FSI™** inserts. These inserts are extremely thin and **Fragile**. During normal use the water may **not mist** like standard inserts, however the cleaning action is just as effective. It is recommended the **Power knob** be set in the Blue area at assure best performance.

Call your **DEALER** or *customer service at Sable Industries Inc.* for further instructions at area **(519) 579-9323**, and a technician will be available to help you with questions, and what to do next.

SCALING TECHNIQUES

Scaling techniques are at the option of the professional using this instrument.

1. Excessive pressure is **not necessary** with **THE AUTOSCALER**. Use a gentle, rapid motion with small overlapping strokes, being certain that the insert is at a *15 degree angle* to the tooth surface.
2. If patient Hypersensitivity is noted, you should *lighten pressure* on the handpiece, *Reduce the Power Level*, or *move from the sensitive tooth* to another, and return to the tooth later. If sensitivity persists, switch to a hand instrument.

WARRANTY INFORMATION

Proof of purchase such as a copy of **ORIGINAL BILL OF SALE** must be submitted to obtain warranty service. Warranty **DOES NOT** include Shipping and Handling charges in either direction. (**See warranty for further details**).

WARRANTY

Sable Industries Inc. (hereafter **SII**) subject to the terms, and conditions herein below expressly set forth in Paragraphs "1" through "4" **HEREBY WARRANTS** that it will repair **THE AUTOSCALER Ultrasonic Scaler** which proves defective by reason of improper workmanship and or materials.

- 1. SIX YEAR:** Parts and Labor coverage on circuit boards, and **TWO YEAR** Parts and Labor coverage on all other components in the device (**with exception of Handpieces**). The longest warranty in the industry.
- 2. ORIGINAL PURCHASER:** This warranty is limited to the **original purchaser** of the unit. When requesting warranty service the following information is needed by **SII**: A copy of the Retail Sales Receipt must be submitted to SII within **2 WEEKS after installation** to validate warranty. A **RETAIL** sales receipt can be used to **VALIDATE** sales/installation date.
- 3. PROPER DELIVERY:** The unit must be shipped, **Freight Prepaid**, or delivered to **SII** to render the services provided hereunder in either its original package, or a similar package affording an equal degree of protection. **Warranty does not include shipping and handling back to the customer.**
- 4. UNAUTHORIZED REPAIR, ABUSE, ETC.:** The unit must not have been previously altered, repaired or serviced by anyone other than a service facility authorized by **SII**. To render such service, the serial number on the unit must not have been altered or removed. The unit must not have been subjected to accident, misuse, abuse, or operated contrary to the operating instructions. **Handpiece assembly** is not warranted if subjected to abuse from excessive pulling/stretching, debris and/or chemical accumulation from water supplies. **Water system** is not warranted by failure due to debris, and/or chemical deposits from water supplies or added medicinal solutions. **Foot switch** is not warranted by failure due to tampering, or operation/storage in a wet environment.

EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, NO OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING WARRANTIES OF MERCHANT ABILITY, AND FITNESS FOR A PARTICULAR PURPOSE, SHALL APPLY TO THIS EQUIPMENT. DAMAGES ARE LIMITED STRICTLY TO REPAIR OR REPLACEMENT OF PARTS. UNDER NO CIRCUMSTANCES SHALL SII BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM THE USE OF THE EQUIPMENT, AND SII NEITHER ASSUMES NOR AUTHORIZES ANY REPRESENTATIVE, OR OTHER PERSON TO ASSUME FOR IT ANY OBLIGATION FOR LIABILITY OTHER THAN SUCH AS IS EXPRESSLY SET FORTH HEREIN.

TABLE OF CONTENTS

3	
30,000HZ (30K) MODEL	5
30K	1, 5, 8, 9
30K STICKER.....	9
A	
AIR SWITCH ACTIVATOR	6
AUTOCLAVE QUESTIONS.....	2
B	
BLUE SHADED AREA	4, 8
C	
CARDIAC PACE MAKER.....	2
CAUTION	2
CAVITRON	2
CORRECT FREQUENCY INSERT	8
CUSTOMER SERVICE	10
D	
DESCRIPTION:.....	3
E	
ELECTRICAL REQUIREMENTS	5
ELECTRONIC GENERATOR/WATER SYSTEM	4
EYE AND FACE PROTECTION.....	3
F	
FEATURES	1
FILTRATION, WATER.....	6
FRAGILE	5
FSI.....	1, 8, 10
FUSE	9
G	
GFI GROUND FAULT INTERRUPT.....	5
GLOVES.....	3
H	
HANDPIECE ASSEMBLY	4
HANDPIECE INSTALLATION	5
HANDPIECE OR INSERT RUNS HOT	9
HANDPIECE PANEL RECEPTACLE.....	4
HANDPIECE SHEATHES	3HIGH POWER
.....	8
HUMMING OR HIGH PITCH SOUND.....	10
I	
INSERT REPLACEMENT	5
INSERTS	1, 4, 5, 8
INSTALLATION / OPERATION.....	6
INTRODUCTION	3
L	
L.E.D.	4
L.E.D. LIGHT	6, 9
N	
NOTE	2, 5, 6, 7
O	
O RING	7
P	
PERIODONTAL INSERTS.....	4, 8
POWER CONTROL KNOB.....	8
POWER CORD.....	6
Q	
QUICK CONNECTOR	3, 5, 6
S	
SCALING TECHNIQUES.....	10
SHUT DOWN PROCEDURE.....	8
SLIMLINE.....	10
SLIMLINE.....	1, 4, 5, 8
SPECIAL PRECAUTIONS.....	2
STERILIZATION	2, 3
STERILIZATION	2, 3, 7
T	
TECHNICAL SERVICE.....	9
TOGGLE SWITCH.....	4
U	
UNIT STOPS WORKING.....	9
W	
WARRANTY	10, 11
WATER	2, 3, 4, 7, 8, 9
WATER KNOB.....	6
WATER KNOB.....	7
WATER KNOB.....	4
WATER, ABSENCE OF	7

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