ULTRASONIC SCALER (25/30K)

P/N 1900005

P/N 1900005(P)

SABLE INDUSTRIES INC.

100 Campbell Avenue - Suite 5, Kitchener, On N2H 4X8 Canada
TEL: 1-800-368-8106 www.sableindustriesinc.com

Table Of Contents

Section Number	Section Title Description of Contents	Page Number
1	Precautions	3
	1.1 Precautions for All Systems	
	1.2 Precautions for Ultrasonic Prophylaxis Procedure	
2	Introduction.....	4
3	Intended Use.....	4
4	Specification	4
5	Contraindications and Warnings.....	5
	5.1 Contraindications	
	5.2 Warnings	
6	Infection Control.....	5
	6.1 General Infection Control Recommendations	
	6.2 Water Supply Recommendations	
7	Installation Instructions	6
	7.1 General Information	
	7.2 Water Line Requirements	
	7.3 Electrical Requirements	
	7.4 Unpacking the System	
	7.5 System Installation	
	7.6 Vertical Installation	
	7.7 Rear Panel Controls / Power Connection	
	7.8 Water Supply Line Connection	
8	Sable Industries Inc® Ultrasonic Scaler Description... ..	7-8
	8.1 System Controls	
	8.2 Handpiece Holder	
	8.3 Handpiece Assembly / Handpiece Sleeve	
	8.4 Foot Switch Controls and Operation	
9	Accessories.....	8
10	Techniques For Use	9
	10.1 Patient Positioning	
	10.2 Performing Ultrasonic Scaling Procedures	
	10.3 Patient Comfort Considerations	
11	System Maintenance And Care.....	10
12	Trouble shooting	11
13	Disposable of Unit	11
14	Disclaimer	11
15	WARRANTY	12
16	Additional Symbols.....	12

Section 1: Precautions

Prior to installation and start-up of the Ultrasonic Scaler, carefully read the instructions provided herein!

1.1 Precautions for All Systems

Do not place the Ultrasonic Scaler on or next to a radiator or other heat source. Excessive heat may damage the Ultrasonic Scaler's electronics. Place the Ultrasonic Scaler where air is free to circulate on all sides and beneath it. Do not cover vents on rear panel.

The Ultrasonic Scaler which can be carried vertical or plate is portable, but must be handled with care when moving.

Equipment flushing and dental water supply system maintenance are strongly recommended. See Section 11: System Maintenance And Care.

Close the water shut-off valve in the dental water supply system every night before leaving the office.

The use of an in-line water filter is recommended.

Never operate the Ultrasonic Scaler without water flowing through the handpiece cable set.

Grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle marked Hospital Only or Hospital Grade.

The insert tip is the applied part intended to contact the patient during use.

1.2 Precautions for Ultrasonic Prophylaxis Procedures

- Like a toothbrush, ultrasonic inserts "wear out" with use. Inserts with just 2 mm of wear lose about 50% of their scaling efficiency. In general, it is recommended that ultrasonic inserts be discarded and replaced after one year of use to maintain optimal efficiency and avoid breakage.
- If excessive wear is noted, or the insert has been bent, reshaped or otherwise damaged, discard the insert immediately.
- Ultrasonic insert tips that have been bent, damaged, or reshaped are susceptible to in-use breakage and should be discarded and replaced immediately.
- Retract the lips, cheeks and tongue to prevent contact with the insert tip whenever it is placed in the patient's mouth.

Section 2: Introduction

Sable Industries Inc®. is an ISO 13485: 2016 certified Company. All UltraSonic Scalers are FDA registered and contain a Health Canada License: 73291.

2.1 Conformance to Standards:

The Sable Industries Inc® Ultrasonic Scaler conforms to Safety Regulations: IEC 60601-1:2005, IEC 60601-1:2005/AMD1:2012, IEC 60601-1:2005/AMD2:2020, EMC Regulation: IEC 60601-1-2:2014 +A1:2020 and Canadian Safety Standards: SPE-3000

2.2 Supplies and Replacement Parts

Contact your local dealer to order supplies or replacement parts. There are no serviceable parts included in this device. Please contact your dealer to acquire all repairing service and technique supports.

CAUTION: EQUIPMENT NOT SUITABLE FOR USE IN THE PRESENCE OF FLAMMABLE ANESTHETIC MIXTURE WITH AIR OR WITH NITROUS OXIDE.

CAUTION: U.S FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A DENTAL PROFESSIONAL.

Section 3: Intended Use

Ultrasonic procedure

The Ultrasonic Scaler is designed for use in prophylaxis treatments.

THIS DEVICE SHOULD BE USED ONLY BY LICENSED DENTAL PROFESSIONALS TRAINED IN THE PROPER USE OF ULTRASONIC SCALING DEVICES.

Section 4: Specifications

Sable Industries Inc® Ultrasonic Scaler:

Multi-voltage unit.

4.1 Electric Voltage:

VOLTAGE:	CURRENT
100VAC-240VAC 50/60Hz	1.2 A (MAX)
WATTAGE CONSUMPTION:	MAX. 85W
WATER PRESSURE:	20-40 psi

Tip vibration frequency: 25/30kHz±2KHz

Half-excursion force: ≤1N

4.3 Operation Environment:

AMBIENT TEMPERATURE	: +10°C ~ +40°C
RELATIVE HUMIDITY	: 30% ~ 75%
ATMOSPHERIC PRESSURE	: 700 hPa ~ 1060hPa

4.2 Dimension:

HEIGHT :	6.5 CM
WIDTH :	23.0 CM
DEPTH :	22.0 CM
WEIGHT :	2.1 KG

4.4 Transportation And Storage Environment:

AMBIENT TEMPERATURE	: - 10°C ~ +70°C
RELATIVE HUMIDITY	: 15%~90%
ATMOSPHERIC PRESSURE	: 500 hPa ~ 1060hPa

Section 5: Contraindications and Warnings

5.1 Contraindications

- Ultrasonic Systems should not be used for restorative dental procedures involving the condensation of amalgam.
- This device is designed to work with all Cavitron® (remark) insert with 25Khz or 30Khz frequency automatically. For optimum performance please use only qualified inserts supplied by Cavitron® or Sable Industries Inc®.
- Do not use this device if the patient or operator is wearing a pacemaker.

Remark: Cavitron® is a registered trademark of Dentsply® International, Inc and is in no way affiliated with Sable Industries Inc®.

5.2 Warnings

- Persons who are fitted with cardiac pacemakers, defibrillators and other active implanted medical devices, have been cautioned that some types of electronic equipment might interfere with the operation of the device. We recommend that the handpiece cable kept at least 6 to 9 inches (15 to 23 cm) away from any device and their leads during use.
- During using the unit, make sure that water is flowing continuously. If the handpiece overheats please check the water supply, and stop using the unit for a while.

Section 6: Infection Control

6.1 General Infection Control

Recommendations

- As with all dental procedures, the use of standard personal protection equipment (i.e., wearing a face mask, eye wear, or face shield, gloves and protective gown) is recommended.
- For maximal operator and patient safety, carefully follow section 11 system maintenance and care information detailed in the operating instruction
- As with high speed handpieces, and other dental devices, the combination of water and ultrasonic vibration from your Ultrasonic Scaler will create aerosols. With proper technique, much of the aerosol dispersion can be effectively controlled and minimized. Please carefully follow the procedural guide lines in this manual regarding the use of your ultrasonic scaler.
- Always flush your Ultrasonic Scaler with highest flow before treatment. Refer to more information in section 10.
- Clean and disinfect the Sterilizable Sleeve of handpiece between patients. The Sterilizable Sleeve of handpiece can be autoclaved up to 135

°C for at least 3 minutes.

Sterilizing:

1. Place handpiece sleeve and pouched ultrasonic insert into a Dynamic Air Removal steam autoclave. After warm-up is completed, operate at a sterilizing temperature and pressure of 273° F/31 psi (134°C/216 kPa) for 4 minutes, followed by a 20-30 minute drying time.
2. To maintain sterility, the insert should remain in the sealed pouch until it is ready for use.

DO NOT USE Cold sterilization solution.

6.2 Water Supply Recommendations

- It is highly recommended that all dental water supply systems should conform to applicable CDC (Centers for Disease Control and Prevention) and ADA (American Dental Association) standards, and that all recommendations be followed in terms of flushing, chemical flushing, and general infection control procedures. See sections 6.1 and 11.

Section 7: Installation Instructions

7.1 General Information

If the installation of your Ultrasonic Scaler is performed by someone other than trained distributor personnel, care should be taken to observe the following requirements and recommendations.

7.2 Water Line Requirements

- The System's water supply line is factory installed. Do not disconnect from the ultrasonic scaler.
- Incoming water supply line pressure to the Ultrasonic Scaler must be 20 psi (172 kPa minimum) to 40 psi (276 kPa) maximum. If your dental water system's supply line pressure is above 60 psi, install a water pressure regulator on the water supply line to your Ultrasonic Scaler.
- A manual shut-off valve on the dental water system supply line should be used so that the water can be completely shut-off when the office is unoccupied.
- A filter in the dental water system supply line is recommended so that any particles in the water supply will be trapped before reaching the Ultrasonic Scaler.
- After the above installations are completed on the dental water supply system, the dental office water line should be thoroughly flushed prior to connection to the System.
- After flushing system verify there are no leaks.

7.3 Electrical Requirements

Refer to Section 4: Specifications.

7.4 Unpacking the System

Carefully unpack your Ultrasonic Scaler and verify that all components and accessories are included:

1. Ultrasonic Scaler main unit with factory installed water supply line and handpiece cable set.
2. Power Cord
3. Foot Pedal.


4. User Manual.

5. Ultrasonic Inserts (Optional).

7.5 System Installation

- The Ultrasonic Scaler Unit is designed for both horizontal and vertical placement on a level surface. Refer to section 7.6 for more information of vertical installation.
- Be sure unit is stable and resting on four feet.
- Placing unit in direct sunlight may discolor plastic housing.

7.6 Vertical Installation

7.7 Rear Panel Controls / Power Connection


- The ON/OFF Control Switch is located on the Front Panel of the Ultrasonic Scaler. (See Section 8.1)
- Verify the ON/OFF Control Switch located on the front panel of the Ultrasonic Scaler is in the OFF position before proceeding.
- Plug the detachable Power Cord into the back of the ultrasonic scaler and into an approved outlet.

7.8 Water Supply Line Connection


Connect the free end of the Ultrasonic Scaler's water supply line to the dental water supply line. Inspect all connections to make certain there are no leaks.

Section 8: Sable Industries Inc® Ultrasonic Scaler Description

8.1 System Controls – Rear Panel



Front Panel

8.2 Handpiece Holder

When your Ultrasonic Scaler is set up horizontal, you can choose to install handpiece holder to the left side of the main unit if you are a left-handed user.



8.3 Sterilizable Sleeve of Handpiece Assembly

The Ultrasonic Scaler is multi-frequency unit compatible with all Cavitron® 25Khz or 30 Khz inserts. The system will automatically detect the insert frequency, no need to switch any button.

For more oral hygiene care, the handpiece sleeve can be dismantled and autoclaved. See below:

Sterilizable Sleeve of Handpiece

Enable to dismantled and autoclave

Fastening Ring

Hold here, rotate
Sterilizable sleeve of
Handpiece to release.

To take off the Sterilizable sleeve of handpiece, please hold the fastening ring, and rotate the sleeve counterclockwise from handpiece cable set assembly. **Note: Avoid rotate or twist the handpiece Cable.**

8.4 Foot Pedal Controls and Operation

The Foot Pedal is a two-position momentary switch, which activates both Ultrasonic energy and water flow at insert tip.

Foot Pedal Released... Both ultrasonic activation and irrigating flow stop.	Foot Pedal Depressed Half Way (1st position) The ultrasonic insert tip is activated and irrigating water flow under normal function.	Foot Pedal fully depressed (2nd position) This activates the "Turbo" mode increasing extra power.

Section 9: Accessories

1. Ultrasonic insert (optional)
2. Power Cord
3. Foot Pedal
4. Handpiece Cable Set

Section 10: Techniques For Use

10.1 Patient Positioning

- For optimal access to both the upper and lower arches, the backrest of the chair should be adjusted to a 45° degree angle. This assures patient comfort and Clinician visibility.
- Have the patient turn his/her head to the right or left. Also position chin up or down depending upon the quadrant and surface being treated. Evacuate irrigate using either a saliva ejector or High Volume Evacuator (HVE).

10.2 Performing Ultrasonic Scaling Procedures

- Note: Refer to the Infection Control Information Booklet supplied with your system for general procedures to be followed at the beginning of each day and between patients.
- The edges of ultrasonic inserts are intentionally rounded so there is little danger of tissue laceration with proper ultrasonic Scaling technique. Whenever the insert tip is placed in the patient's Mouth, the lips, cheek and tongue should be retracted to prevent Accidental prolonged contact with the activated tip.
- Hold the empty handpiece in an upright position. Activate the Foot Control until fluid exits.
- Lubricate the rubber O-ring on the insert with water before placing it into the handpiece. Fully seat insert with a gentle push-twist motion. DO NOT FORCE IT INTO PLACE.
- Activate the System. Hold the handpiece over a sink or drain. Check spray temperature to verify fluid is reaching the working end of the insert tip. Adjust the water cooler irrigate. Control knob to ensure adequate flow for the selected Power setting. Greater flow settings provide cooler irrigation.
- It may be necessary to adjust water flow larger under "Turbo" mode (Foot pedal fully depressed) so adequate fluid will be available to cool tip and tooth interface.

- In general, it is suggested a "feather-light-touch" be used both supra and subgingivally. The motion of the activated tip and acoustic effects of the irrigating fluid, in most cases, is adequate to remove even the most tenacious calculus.
- Periodically check the ultrasonic insert for wear with the Insert Efficiency Indicator.
- The use of a saliva ejector or High Volume Evacuator (HVE) is recommended during all procedures.
- Set the System's Power Intensity display to the lowest power setting for the application and the selected insert.
- If water leakage found in handpiece, replace sleeve or o-ring on handpiece cable set assembly to eliminate.

Note: If the Power Intensity control is set to minimal, it might get chance you can not activate the system well, like no water flow or not vibrating. In that case, please set the Power Intensity control up to like level 3, or 4, then activate the Foot Control and then turn the Power Intensity Adjustment Knob to your preferred setting.

10.3 Patient Comfort Considerations

Reasons for sensitivity

- Incorrect tip placement. Point should be directed away from root surfaces.
- Not keeping tip in motion on tooth. Do not allow the insert to remain in a static position on any one area of the tooth. Change the insert's path of motion.
- Applying pressure. Use extremely light grasp and pressure, especially on exposed cementum.
- If sensitivity persists, decrease power setting and/or move from the sensitive tooth to another and then return.

Section 11: System Maintenance and Care

Daily Maintenance

It is recommended that you perform the following maintenance procedures to help minimize bio-film formation in the water path of your Ultrasonic Scaler which could affect the water flow to the ultrasonic insert, and scaling performance.

Start-Up Procedures at the beginning of the day:

1. Open the manual shut-off valve on the dental office water supply system.
2. Turn the System ON using the Power ON/OFF switch. (see illustration on page 7) Verify the Intensity display is on.
3. Set the Intensity display to minimum setting.
4. Set the Water Control knob to maximum.
5. Hold the Sterilizable sleeve of Handpiece (without an insert installed) upright over a sink or drain. Activate the Foot Control and flush the water line for at least 2 minutes.
6. Place a sterilized insert into the sterilizable sleeve of Handpiece and set the water control knob to your preferred operating position.

Between Patients:

1. Remove ultrasonic insert and sterilizable sleeve of handpiece used, clean and sterilize.
2. Clean and disinfect the surfaces of the cabinet, Power Cord, Handpiece cable set assembly*, control knobs, Foot pedal and cable assembly by applying an approved non-immersion type disinfectant solution* carefully following the instructions provided by the disinfectant solution manufacturer. To clean system, generously spray disinfectant solution on a clean towel and wipe all surfaces. Discard used towel. To disinfect system, generously spray disinfectant on a clean towel and wipe all surfaces. Allow disinfectant solution to air dry. Do not spray disinfectant solution directly on the Ultrasonic Scaler.
3. Place a sterilized handpiece sleeve. Set power to minimum. Hold the handpiece over a sink or drain

and flush the water line at maximum water flow for 30 seconds.

4. When ready, place a sterilized insert into the handpiece.
5. Please consider use a FDA approved sheath to for the cable and entire sterilizable sleeve of handpiece or at least to cover the sterilizable handpiece from the fastening ring to cable when sterilizable sleeve of handpiece is sterilized between patients.

Shut-Down Procedures at the end of the day:

1. Remove ultrasonic insert and handpiece sleeve used, clean and sterilize.
2. Turn the System OFF.
3. Clean and disinfect the surfaces of the cabinet, Power Cord, Handpiece cable set assembly, control knobs, Foot pedal cable assembly by applying an approved non-immersion type disinfectant solution* carefully following the instructions provided by the disinfectant solution manufacturer. To clean system, generously spray disinfectant solution on a clean towel and wipe all surfaces. Discard used towel. To disinfect system, generously spray disinfectant on a clean towel and wipe all surfaces. Allow disinfectant solution to air dry. Do not spray disinfectant solution directly on the ultrasonic scaler.
4. Close the manual shut-off valve on the dental water supply system.

***Note:** Those Water-based disinfection solutions are preferred, due some alcohol-based disinfectant solutions may be harmful and may discolor plastic materials.

***Note:** The electric wire winding covered by heat shrinkage tube which protected by sterilizable sleeve of handpiece is very sensitive to the disinfectant solution and water. After the disinfection, wipe surface of shrinkage tube with a slightly damp cloth and dry thoroughly before use.

Section 12: Trouble shooting

Although service and repair of the Ultrasonic Scaler should be performed by dealer personnel, the following are some basic trouble shooting procedures that will help avoid unnecessary service calls. Generally, check all lines and connections to and from the System, a loose plug or connection will often create problems. Check the settings on the System's knobs.

12.1 Troubleshooting Guide

Problem: System will not operate: (Intensity display is not lit.)

1. Check that the Power switch is in the ON position, and that the detachable Power Cord is fully seated in the receptacle on back of System.
2. Check that the System's three-prong plug is fully seated in an appropriate AC receptacle, and that AC current is present.

(Intensity display is lit.)

1. Check that the Foot Pedal Connector is fully seated in the Foot Pedal Receptacle on the back of the System.

System operates:

(No water flow to insert.)

1. Assure that water control is properly adjusted.
2. Check that water supply control valve(s) (dental office water supply) are open.

(Insert stops vibrating)

1. Deactivate foot pedal.
2. Verify insert is in good condition.
3. Depress foot pedal to try again.

Section 13: Disposal of Unit

- Keep original packaging until the Ultrasonic Scaler is to be disposed of permanently. You can use it for shipping or storing your Ultrasonic Scaler at any time.
- Dispose of the Ultrasonic Scaler in accordance with local and national laws.

Section 14: Disclaimer

SABLE considers itself responsible for the effects on safety, reliability and performance of this product only if:

- Assembly operations, extensions, re-adjustments, modifications or repairs are carried out by persons authorized by SABLE.
- The electrical installation of the relevant room complies with the requirements.
- The equipment is used in accordance with these instructions for use.

Section 15: Warranty

15.1 Malfunction

SABLE hereby warrants that for a period of two years from the delivery date, this device shall be free from defects in material and workmanship. In case the machine is found malfunctioned under normal use, SABLE will offer service of free maintenance and parts for replacement.

15.2 Repair

Repairs must be only carried out by an authorized SABLE engineer/dealer. If repairs during warranty period are not carried out by an authorized engineer/dealer, warranty will expire immediately.

15.3 Warranty Exception

The warranty stated herein is the sole warranty applicable to SABLE products. SABLE expressly disclaims the liability for warranty even within warranty period, if

- (1) Damages caused by natural disaster.
- (2) Operator's fault or wrong operation.
- (3) Application use other than curing light-cured material purpose.
- (4) A malfunction or damage caused by repair, adjustment, modification which is not carried out by Sable authorized technicians/dealers.
- (5) A malfunction caused by abnormal power source or voltage.
- (6) It is a consumption part.

Section 16: Additional Symbols

BF Type.

Attention, refer to accompanying documents.

Power off

Power on

ANNEX I

Guidance and manufacturer's declaration-electromagnetic emissions		
The SABLE Ultrasonic Scaler is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The Ultrasonic Scaler uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Compliance	The Ultrasonic Scaler is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration-electromagnetic immunity

The Ultrasonic Scaler is intended for use in the electromagnetic environment specified below.

The customer or the user of the Ultrasonic Scaler should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD) IEC 61000-4-2	±8kV contact ±15kV air	±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output	±2kV for power supply lines ±1kV for input/output	Mains power quality should be that of a typical commercial or hospital environment.
Surge immunity IEC 61000-4-5	±1kV line(s) to line(s) ±2kV line(s) to earth	±1kV line(s) to line(s) ±2kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT(>95% dip in UT) for 0.5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles <5% UT(>95% dip in UT) for 5secs	<5% UT(>95% dip in UT) for 0.5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles <5% UT(>95% dip in UT) for 5secs	Mains power quality should be that of a typical commercial or hospital environment. If the user requires continued operation during power mains interruptions, it is recommended that the scaler be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz)magnetic field EN 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration-electromagnetic immunity			
The Ultrasonic Scaler is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultrasonic Scaler should assure that is used in such and environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC61000-4-6	3Vrms 150 kHz to 80MHz	3Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the MS-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80MHz to 800MHz $d = 2.3\sqrt{P}$ 800MHz to 2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC61000-4-3	3V/m 80MHz to 2.5 GHz	3V/m	
Electrical fast transient/burst EN 61000-4-4	±2kV for power supply lines ±1kV for input/output	±2kV for power supply lines ±1kV for input/output	Mains power quality should be that of a typical commercial or hospital environment.
NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Ultrasonic Scaler is used exceeds the applicable RF compliance level above, the Ultrasonic Scaler should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Ultrasonic Scaler.			
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distance between portable and mobile RF communications equipment and the Ultrasonic Scaler

The Ultrasonic Scaler is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Ultrasonic Scaler can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Ultrasonic Scaler as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d=1.2\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2.5GHz $d=2.3\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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: snossəp

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The diagram shows the assembly of the handle and the torque joint. The handle is labeled "pièce à main" (handle piece) and "Assemblage de l'ensemble de main" (assembly of the handle set). The torque joint is labeled "Joint torique" (torque joint). The handle is shown being inserted into the torque joint. The text "permet de démonter et l'autoclaver" (allows for disassembly and autoclaving) is written vertically on the left side of the diagram.

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