

Dental Ops	Restore		Daily Maintenance <i>*These instructions based on a 4-day work week*</i>		Twice Weekly Maintenance	
	Scoops / Grams	Water - oz. & mL	Scoops / Grams	Water - oz. & mL	Scoops / Grams	Water - oz. & mL
1	1/2 Scoop (2.5g)	16oz. / 473mL	1/2 Scoop (2.5g)	16oz. / 473mL	1 Scoop (5g)	16oz. / 473mL
2	1 Scoop (5g)	32oz. / 946mL	1 Scoop (5g)	32oz. / 946mL	2 Scoops (10g)	32oz. / 946mL
3	1 1/2 Scoops (7.5g)	48oz. / 1.419L	1 1/2 Scoops (7.5g)	48oz. / 1.419L	3 Scoops (15g)	48oz. / 1.419L
4	2 Scoops (10g)	64oz. / 1.892L	2 Scoops (10g)	64oz. / 1.892L	4 Scoops (20g)	64oz. / 1.892L
5	2 1/2 Scoops (12.5g)	80oz. / 2.365L	2 1/2 Scoops (12.5g)	80oz. / 2.365L	5 Scoops (25g)	80oz. / 2.365L
6	3 Scoops (15g)	96oz. / 2.838L	3 Scoops (15g)	96oz. / 2.838L	6 Scoops (30g)	96oz. / 2.838L
7	3 1/2 Scoops (17.5g)	112oz. / 3.311L	3 1/2 Scoops (17.5g)	112oz. / 3.311L	7 Scoops (35g)	112oz. / 3.311L
8	4 Scoops (20g)	128oz. / 3.784L	4 Scoops (20g)	128oz. / 3.784L	8 Scoops (40g)	128oz. / 3.784L
9	4 1/2 Scoops (22.5g)	144oz. / 4.257L	4 1/2 Scoops (22.5g)	144oz. / 4.257L	9 Scoops (45g)	144oz. / 4.257L
10	5 Scoops (25g)	160oz. / 4.730L	5 Scoops (25g)	160oz. / 4.730L	10 Scoops (50g)	160oz. / 4.730L
11	5 1/2 Scoops (27.5g)	176oz. / 5.203L	5 1/2 Scoops (27.5g)	176oz. / 5.203L	11 Scoops (55g)	176oz. / 5.203L
12	6 Scoops (30g)	192oz. / 5.676L	6 Scoops (30g)	192oz. / 5.676L	12 Scoops (60g)	192oz. / 5.676L
13	6 1/2 Scoops (32.5g)	208oz. / 6.149L	6 1/2 Scoops (32.5g)	208oz. / 6.149L	13 Scoops (65g)	208oz. / 6.149L
14	7 Scoops (35g)	224oz. / 6.622L	7 Scoops (35g)	224oz. / 6.622L	14 Scoops (75g)	224oz. / 6.622L
15	7 1/2 Scoops (37.5g)	240oz. / 7.095L	7 1/2 Scoops (37.5g)	240oz. / 7.095L	15 Scoops (75g)	240oz. / 7.095L
16	8 Scoops (40g)	256oz. / 7.568L	8 Scoops (40g)	256oz. / 7.568L	16 Scoops (80g)	256oz. / 7.568L
17	8 1/2 Scoops (42.5g)	272oz. / 8.041L	8 1/2 Scoops (42.5g)	272oz. / 8.041L	17 Scoops (85g)	272oz. / 8.041L
18	9 Scoops (45g)	288oz. / 8.514L	9 Scoops (45g)	288oz. / 8.514L	18 Scoops (95g)	288oz. / 8.514L
19	9 1/2 Scoops (47.5g)	304oz. / 8.987L	9 1/2 Scoops (47.5g)	304oz. / 8.987L	19 Scoops (95g)	304oz. / 8.987L
20	10 Scoops (50g)	320oz. / 9.460L	10 Scoops (50g)	320oz. / 9.460L	20 Scoops (100g)	320oz. / 9.460L

Bio-Pure should be mixed with warm water and used at the end of the work day

