


NEW! Autoclave Gaskets



1700301
SciCan Statim 2000
Cassette Type Gasket Seal



1700302
SciCan Statim 5000
Cassette Type Gasket Seal



1700307
Midmark M11 Door Type
Gasket Seal with Dam

3 Wire Buck-Boost Transformer for Easy Set up

Critical Dental Equipment like Compressors, Vacuums, Sterilizers, Pan-Ceph's, and Intraoral X-Rays require a stable supply voltage to maintain optimum performance and prevent premature failure.



Sable Buck-Boost Transformers features:

- ✓ Steps up or steps down incoming voltage.
- ✓ Made in Canada. UL & CSA approved.
- ✓ Front accessible separate high and low voltage lead wires.
- ✓ Valid for single and three phase -
(for three phase 2 Buck-Boost are required).
- ✓ 10 Year Warranty.




3 Wire Set Up Buck-Boost - 3 Sizes to Choose From:

- 3001015 5.0KVA – Single Head compressors or Suction units (up to 15 amp and under).
- 3001025 7.5KVA – Dual Head Compressors or Dual Head Suction units (up to 20 amp and under).

Also available: Original 8 Wire Buck-Boost:

- 3001000 3.5KVA – Small applications such as sterilizers and other assorted table top items (up to 10 amp and under)
- 3001010 5.0KVA – Single Head compressors or Suction units (up to 15 amp and under).
- 3001020 7.5KVA – Dual Head Compressors or Dual Head Suction units (up to 20 amp and under).

Silicone Tubing



High quality silicone Euro type tubings complete with 10-32 to metric barbs for installation in North American dental units. 2100401, 2100402 comes with Quick Disconnect

- 2100301 5 FT, 4-Hole Gray Silicone Tubing
- 2100302 5 FT, 6-Hole Gray Silicone Tubing
- 2100304 7 FT, 4-Hole Gray Silicone Tubing
- 2100305 7 FT, 6-Hole Gray Silicone Tubing
- 2100401 4-Hole A-dec* Type Gray Silicone Tubing
- 2100402 6-Hole A-dec* Type Gray Silicone Tubing



6 PIN



5 Hole

6 Hole, Fiber Optic, Straight Tubing

- | | |
|--------------------------|---------------------------|
| 2100103 7 Foot Gray | 2100104 7 Foot Surf |
| 2100105 10 Foot Gray | 2100106 10 Foot Surf |
| 2100108 7 Foot Sterling | 2100109 10 Foot Sterling |
| 2100110 7 Foot Dark Surf | 2100111 10 Foot Dark Surf |

5 Hole, Fiber Optic, Straight Tubing with Bulb

- 2100201 7 Foot Sterling
- 2100202 7 Foot Gray

Replacement Bulb for all Sable 5-Hole Tubing, please order: 1000815



Light Board

- Lamp/LED Light Control Board powers all makes including STAR* bulbs
- Power Transformer 9VDC 2.2Amps

Light Control Kits:

- 2300001 6 hole, 7ft. - Gray** Tubing
- 2300002 6 hole, 10ft. - Gray** Tubing
- 2300003 5 hole, 7ft. - Gray** Tubing w/Bulb

Ask us about more color options