


# ***SABLE Rotamax***

## ***45° High Speed Handpieces***

- Torque Head
- 20 Watts Power
- 350,000 RPMs
- Rear Venting Back Cap
- Ceramic German Bearings
- Titanium Plated Finish
- Fiber Optic with a Kavo® Type Connection
- Midwest® Type Fixed 4 hole Non-Fiber Optic Connection
- 1 Year Warranty



2000023 Rotamax™ PRO 45KL - Fiber Optic

2000024 Rotamax™ PRO 45 - Non-Optic

The following are registered trademarks of the indicated companies: Kavo - Kaltenbach & Voigt GmbH & Co., DENTSPLY International Inc (Midwest). Sable Industries Inc. is in no way affiliated with the above companies.

CELEBRATING  
**25**  
YEARS OF  
EXCELLENCE

**SABLE**  
INDUSTRIES INC.

[WWW.SABLEINDUSTRIESINC.COM](http://WWW.SABLEINDUSTRIESINC.COM)

[f](#) [@](#) [in](#) 1-800-368-8106